




**2mg rProtein Guaranteed, Mammalian Cell Expression,
21 to 30 days, For USD 2500 only!**

This promotion includes:

- cDNA cloning into mammalian expression vectors
- Plasmid DNA transfection and protein expression in HEK293 cells
- mAb purification using His-tag or protein A/G affinity purification
- Risk free after feasibility evaluation
- Scale up bulk production for **BIG discount**

Within 30days, more than 2 mg of purified protein(>90% purity by SDS-PAGE) and a comprehensive project report will be delivered to customer, two or three weeks needed for gene synthesizing with extra cost if gene sequence information available only.

Process flow of HEK Max[®] transient expression production



For more information please send email to our experts via:
project@acrobiosystems.com(preferable)

Or contact us directly by:
Tel: +1 301-825-5518 (North America)
+44 20-8123-6518 (EU)
+ 86 400-685-2052 (Asia and Pacific)
Fax : +1 301-500-2214