



Product Details

Human Laminin 511 Protein (EasyScale) is a recombinant human laminin 511 E8 (LN511-E8) protein that provides a defined surface for in vitro feeder-free culture of multiple human pluripotent stem cells (PSCs).

As published by Takamichi Miyazaki et al., the LN511-E8 variant of laminin 511 shows higher efficiency for supporting the adhesion of dissociated cells than did wild-type laminin 511 which makes it a cost-effective choice.

LN511-E8 has been proven to maintain normal growth characteristics and stemness in multiple PSC lines without simultaneous differentiation, which includes ESC, iPSC, MSC, etc. In addition, LN511-E8 has been demonstrated to support PSC growth for > 10 passages without any signs of karyotypic abnormalities and to maintain the ability of PSCs to differentiate into all three germ line lineages.

Importantly, our innovative laminin 511 protein (EasyScale) could meanwhile support coating method and uncoating method (directly mixed the protein with cell suspensions without coating plates).

Suitable for automated systems

Easy to use and without coating plates to reduce steps.

Cost effective

Half of amount ($0.25 \mu\text{g}/\text{cm}^2$) used in coating methods and save money.

Stemness maintenance

No spontaneous differentiation is observed after several passages of hPSC culture.

Lot-to-lot consistency

Produced from a stable cell line, robust purification process, stringent QC.

Ready to scale supply

A new cGMP-compliant facility under construction.

Key parameter

Purity (SDS PAGE)	> 95%
Mycoplasma Test	Negative
Sterility Test	Negative
Endotoxin Test	< 0.01 EU per μg

Formulation

Supplied as $0.2 \mu\text{m}$ filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Shipping

This product is supplied and shipped with dry ice, please inquire the shipping cost.

Storage

For long term storage, the product should be stored at liquid state at -70°C .


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- $2-8^\circ\text{C}$ for 24 months under sterile condition;
- -70°C for 36 months.

Bioactivity-Stem Cell Culture

[View Protocol](#)



Robust, homologous, monolayer expansion of single cell hPSCs on Human Laminin 511 Protein, premium grade (EasyScale) (Cat. No. LA1-H5213).

Discounts, Gifts,
and more!





Human Laminin 511 Protein, premium grade (EasyScale)

Catalog # LA1-H5213




BIOSYSTEMS
Acro



hPSCs cultured on Human Laminin 511 Protein, premium grade (EasyScale) (Cat. No. LA1-H5213) maintain pluripotency, exhibiting high expression of OCT4, SOX2, and SSEA4 as confirmed by FACS and immunofluorescent staining.

Bioactivity-SPR



Human Laminin 511 Protein, premium grade (EasyScale) (Cat. No. LA1-H5213) immobilized on CM5 Chip can bind Human ITGA6&ITGB1 Heterodimer Protein, His Tag&Tag Free (Cat. No. IT1-H52W7) with an affinity constant of 9.98 nM as determined in a SPR assay (Biacore 8K) (QC tested).

Clinical and Translational Updates

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com